

FUCA1 antibody

Catalog Number: 110790

Product name

FUCA1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FUCA1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human heart tissue, human kidney tissue, human liver cancer tissue, human placenta tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of FUCA1 recombinant protein (Accession Number: BC017338). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human liver using Catalog No:110790(FUCA1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:110790(FUCA1 antibody) at dilution of 1:50 (under 40x lens)