

**Product name**

GALNT7 antibody

Polyclonal

**Ig Type**

Rabbit IgG

**Specificity**

Human, Mouse, Rat; other species not tested.

**Applications**

ELISA, WB

**Antibody description**

GALNT7 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, human brain tissue. Observed molecular weight by Western-blot: 75kd


**Dilutions**

Recommended Dilution:

WB: 1:500-1:5000

**Preparation**

This antibody was obtained by immunization of GALNT7 recombinant protein (Accession Number: NM\_017423). Purification method: Antigen affinity purified.

**Validations****Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110848(GALNT7 antibody) at dilution of 1:800

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**