

GCK antibody

Catalog Number: 110905

Product name

GCK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GCK Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, mouse pancreas tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 50-52kd

Preparation

This antibody was obtained by immunization of GCK recombinant protein (Accession Number: NM_000162). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110905(GCK antibody) at dilution of 1:500



IP Result of anti-GCK (IP:Catalog No:110905, 4ug; Detection:Catalog No:110905 1:800) with mouse liver tissue lysate 4000ug.