

GEMIN4 antibody



Catalog Number: 110938

Product name

GEMIN4 antibody

Specificity

Human; other species not tested.

Antibody description

GEMIN4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Positive IP detected in HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 120-130 kDa

Preparation

This antibody was obtained by immunization of GEMIN4 recombinant protein (Accession Number: NM_015721). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110938(GEMIN4 antibody) at dilution of 1:100



IP result of anti-GEMIN4 (Catalog No:110938 for IP and Detection) with 293 cell lysate.