

HAAO antibody

Catalog Number: 111328

Product name

HAAO antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HAAO Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, mouse liver tissue, rat liver tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of HAAO recombinant protein (Accession Number: NM_012205). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111328(HAAO antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:111328(HAAO antibody) at dilution of 1:100 (under 10x lens)