

HPS1 antibody



Catalog Number: 111451

Product name

HPS1 antibody

Specificity

Human; other species not tested.

Antibody description

HPS1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HEK-293 cells, K-562 cells, Raji cells, Transfected HEK-293 cells, Y79 cells. Observed molecular weight by Western-blot: 70-79 kDa

Preparation

This antibody was obtained by immunization of HPS1 recombinant protein (Accession Number: NM_182639). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111451(HPS1 antibody) at dilution of 1:800