

HNRNPA0 antibody



Catalog Number: 111499

Product name

HNRNPA0 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

HNRNPA0 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, fetal human brain tissue, human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 32 kDa, 34 kDa

Preparation

This antibody was obtained by immunization of HNRNPA0 recombinant protein (Accession Number: NM_006805). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog

No:111499(HNRNPA0 Antibody) at dilution of 1:600