

# HNRNPA0 antibody



Catalog Number: 111499

---

## Product name

HNRNPA0 antibody

## Immunogen

## Specificity

Human, Mouse; other species not tested.

## Antibody description

HNRNPA0 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, fetal human brain tissue, human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 32 kDa, 34 kDa

## Preparation

This antibody was obtained by immunization of HNRNPA0 recombinant protein (Accession Number: NM\_006805). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111499(HNRNPA0 Antibody) at dilution of 1:600