

IL 1 alpha antibody



Catalog Number: 111652

Product name

IL 1 alpha antibody

Immunogen

Specificity

Human, Mouse; other species not tested.

Antibody description

IL 1 alpha Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human placenta tissue, Jurkat cells, mouse lung tissue. Observed molecular weight by Western-blot: 31 kDa, 17 kDa

Preparation

This antibody was obtained by immunization of IL 1 alpha recombinant protein (Accession Number: NM_000575). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111652(IL1A Antibody) at dilution of 1:600