

# HTRA4 antibody



Catalog Number: 111716

## Product name

HTRA4 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

human, mouse, rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

HTRA4 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 51kd

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of HTRA4 recombinant protein (Accession Number: NM\_153692). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:111716(HTRA4 antibody) at dilution of 1:500

## Clonality