

HUWE1 antibody



Catalog Number: 111729

Product name

HUWE1 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HUWE1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293. Observed molecular weight by Western-blot: 482 kDa

Preparation

This antibody was obtained by immunization of HUWE1 recombinant protein (Accession Number: BC002602). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: N/A

Validations



WB result of HUWE1 antibody (Catalog No:111729, 1:1000) with 293 cells by Dr. Emilia Galperin, University of Kentucky.