

IL-36 γ / IL-1F9 antibody



Catalog Number: 111742

Product name

IL-36 γ / IL-1F9 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200


Immunogen

[Human IL1F9 Recombinant protein \(GST tag\)](#)

Validations

Specificity

Human; other species not tested.




Antibody description

IL-36 γ / IL-1F9 Rabbit Polyclonal antibody.
Positive IHC detected in human oesophagus tissue, human colon cancer tissue. Positive WB detected in Transfected HEK-293 cells.

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111742(IL1F9 Antibody) at dilution of 1:500

Preparation

This antibody was obtained by immunization of IL-36 γ / IL-1F9 recombinant protein (Accession Number: NM_019618). Purification method: Antigen affinity purified.



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:111742(IL1F9 Antibody) at dilution of 1:50 (under 10x lens)

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:111742(IL1F9 Antibody) at dilution of 1:50

IL-36 γ / IL-1F9 antibody



Catalog Number: 111742

(under 40x lens)