

# IL-36γ / IL-1F9 antibody



Catalog Number: 111742

## Product name

IL-36γ / IL-1F9 antibody

## Specificity

Human; other species not tested.

## Antibody description

IL-36γ / IL-1F9 Rabbit Polyclonal antibody.  
Positive IHC detected in human oesophagus tissue, human colon cancer tissue. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of IL-36γ / IL-1F9 recombinant protein (Accession Number: NM\_019618). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111742(IL1F9 Antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:111742(IL1F9 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:111742(IL1F9 Antibody) at dilution of 1:50 (under 40x lens)