

# IL9 antibody



Catalog Number: 111750

## Product name

IL9 antibody

IHC: 1:20-1:200


## Specificity

Human; other species not tested.

## Validations

## Antibody description

IL9 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human lung tissue, human ovary tissue, human skin tissue. Positive WB detected in U-937 cells. Observed molecular weight by Western-blot: 16 kDa




## Preparation

This antibody was obtained by immunization of IL9 recombinant protein (Accession Number: NM\_000590). Purification method: Antigen affinity purified.

U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:111750(IL9 antibody) at dilution of 1:300

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human kidney using Catalog No:111750(IL9 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human kidney using Catalog No:111750(IL9 antibody) at dilution of 1:50 (under 40x lens)