

IL9 antibody

Catalog Number: 111750

Product name

IL9 antibody

IHC: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

IL9 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human lung tissue, human ovary tissue, human skin tissue. Positive WB detected in U-937 cells. Observed molecular weight by Western-blot: 16 kDa




Preparation

This antibody was obtained by immunization of IL9 recombinant protein (Accession Number: NM_000590). Purification method: Antigen affinity purified.

U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:111750(IL9 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.




Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human kidney using Catalog No:111750(IL9 antibody) at dilution of 1:50 (under 10x lens)

Clonality

Polyclonal



Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human kidney using Catalog No:111750(IL9 antibody) at dilution of 1:50 (under 40x lens)