

# Interferon alpha antibody



Catalog Number: 111868

## Product name

Interferon alpha antibody

## Specificity

Human; other species not tested.

## Antibody description

Interferon alpha Rabbit Polyclonal antibody.  
Positive WB detected in Recombinant protein.

## Preparation

This antibody was obtained by immunization of Interferon alpha recombinant protein (Accession Number: NM\_024013). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:111868(IFNA1 Antibody) at dilution of 1:3000