

ITIH3 antibody



Catalog Number: 111948

Product name

ITIH3 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ITIH3 Rabbit Polyclonal antibody. Positive IHC detected in human stomach cancer tissue. Positive IP detected in human plasma tissue. Positive WB detected in human serum tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of ITIH3 recombinant protein (Accession Number: NM_002217). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human serum tissue were subjected to SDS PAGE followed by western blot with Catalog No:111948(ITIH3 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human stomach cancer slide using Catalog No:111948(ITIH3 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human stomach cancer slide using Catalog

ITIH3 antibody



Catalog Number: 111948

No:111948(ITIH3 Antibody) at dilution of 1:50

IP Result of anti-ITIH3 (IP:Catalog No:111948, 4ug; Detection:Catalog No:111948 1:1000) with human plasma tissue lysate 4000ug.

