

ITK antibody



Catalog Number: 111950

Product name

ITK antibody

Immunogen

[Human ITK Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

ITK Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse spleen tissue. Observed molecular weight by Western-blot: 68-72kd

Preparation

This antibody was obtained by immunization of ITK recombinant protein (Accession Number: NM_005546). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111950(ITK antibody) at dilution of 1:300