

KIDINS220 antibody



Catalog Number: 111998

Product name

KIDINS220 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KIDINS220 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, human brain tissue, Raji cells. Observed molecular weight by Western-blot: 186-197 kDa

Preparation

This antibody was obtained by immunization of KIDINS220 recombinant protein (Accession Number: NM_001348736). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111998(KIDINS220 antibody) at dilution of 1:1000