

KGA/GAC antibody



Catalog Number: 112033

Product name

KGA/GAC antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

KGA/GAC Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse brain tissue, mouse kidney tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human gliomas tissue, human brain tissue, human kidney tissue. Observed molecular weight by Western-blot: 65kd, 58kd

Preparation

This antibody was obtained by immunization of KGA/GAC recombinant protein (Accession Number: NM_014905). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, IP, WB


Dilutions

Recommended Dilution:


Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112033 (GLS antibody) at dilution of 1:100 (under 10x lens)



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112033 (GLS Antibody) at dilution of 1:600



KGA/GAC antibody



Catalog Number: 112033

IP Result of anti-GLS (IP:Catalog No:112033, 4ug; Detection:Catalog No:112033 1:500) with mouse brain tissue lysate 8000ug.