

# LAYN antibody



Catalog Number: 112159

## Product name

LAYN antibody

## Immunogen

[Human LAYN Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

LAYN Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue, MCF7 cells. Observed molecular weight by Western-blot: 43kd

## Preparation

This antibody was obtained by immunization of LAYN recombinant protein (Accession Number: NM\_178834). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112159(LAYN antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using LAYN antibody Catalog No:112159 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).