

LGALS3BP antibody



Catalog Number: 112205

Product name

LGALS3BP antibody

Immunogen

Specificity

Human; other species not tested.

Antibody description

LGALS3BP Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, HEK-293 cells, HeLa cells. Positive IP detected in HepG2 cells, HEK-293 cells. Observed molecular weight by Western-blot: 65-90 kDa

Preparation

This antibody was obtained by immunization of LGALS3BP recombinant protein (Accession Number: NM_005567). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:112205(LGALS3BP antibody) at dilution of



1:1000

IP Result of anti-LGALS3BP (IP:Catalog No:112205, 4ug; Detection:Catalog No:112205 1:300) with HepG2 cells lysate 3480ug.