

# LGI2 antibody

Catalog Number: 112207

---

## Product name

LGI2 antibody

## Immunogen

[Human LGI2 Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

LGI2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 62kd

## Preparation

This antibody was obtained by immunization of LGI2 recombinant protein (Accession Number: NM\_018176). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112207(LGI2 antibody) at dilution of 1:800