

# Lumican antibody



Catalog Number: 112394

## Product name

Polyclonal

Lumican antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

Lumican Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 65-90 kDa

## Dilutions


Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of Lumican recombinant protein (Accession Number: NM\_002345). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112394(LUM antibody) at dilution of 1:800

## Clonality