

MBP-Tag antibody

Catalog Number: 112509

Product name

MBP-Tag antibody

Specificity

Recombinant protein; other species not tested.

Antibody description

MBP-Tag Rabbit Polyclonal antibody. Positive WB detected in lane1: MBP-NOL6 64kd; lane2 MBP 42kd, Recombinant protein.

Preparation

This antibody was obtained by immunization of MBP-Tag recombinant protein. Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



lane1: MBP-NOL6 64kd; lane2 MBP 42kd were subjected to SDS PAGE followed by western blot with Catalog No:117331(MBP-Tag Antibody) at dilution of 1:3000