

# MFAP4 antibody



Catalog Number: 112525

## Product name

MFAP4 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MFAP4 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tissue, human lung tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Positive WB detected in human kidney tissue, mouse heart tissue, mouse lung tissue, rat lung tissue. Observed molecular weight by Western-blot: 36kd

## Preparation

This antibody was obtained by immunization of MFAP4 recombinant protein (Accession Number: NM\_002404). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


## Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:112525(MFAP4 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human ovary using Catalog No:112525(MFAP4 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human ovary using Catalog No:112525(MFAP4

# MFAP4 antibody



Catalog Number: 112525

---

antibody) at dilution of 1:100 (under 40x lens)