

MCRS1 antibody



Catalog Number: 112561

Product name

MCRS1 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MCRS1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, human liver tissue, human spleen tissue, Jurkat cells, mouse spleen tissue, mouse uterus tissue, rat spleen tissue. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of MCRS1 recombinant protein (Accession Number: NM_006337). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112561(MCRS1 antibody) at dilution of 1:500