

MEA1 antibody



Catalog Number: 112576

Product name

MEA1 antibody

Immunogen

[Human MEA1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEA1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in Jurkat cells, mouse testis tissue. Observed molecular weight by Western-blot: 34kd and 36kd

Preparation

This antibody was obtained by immunization of MEA1 recombinant protein (Accession Number: NM_014623). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112576(MEA1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:112576(MEA1 antibody) at dilution of 1:50 (under 40x lens)