

MIG antibody



Catalog Number: 112620

Product name

Polyclonal

MIG antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

MIG Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells, human stomach tissue.

Dilutions


Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of MIG recombinant protein (Accession Number: NM_002416). Purification method: Antigen Affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112620(CXCL9 Antibody) at dilution of 1:1000

Clonality