

# MEOX2 antibody



Catalog Number: 112654

## Product name

MEOX2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MEOX2 Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of MEOX2 recombinant protein (Accession Number: NM\_005924). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:112654(MEOX2 Antibody) at dilution of 1:600



IP Result of anti-MEOX2 (IP:Catalog No:112654, 4ug; Detection:Catalog No:112654 1:300) with human placenta tissue lysate 1520ug.