

MT3 antibody

Catalog Number: 112762

Product name

MT3 antibody

Immunogen

[Human MT3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MT3 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human kidney tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of MT3 recombinant protein (Accession Number: NM_005954). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:
IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human brain using Catalog No:112762(MT3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:112762(MT3 antibody) at dilution of 1:50 (under 40x lens)