

MTP18 antibody



Catalog Number: 112793

Product name

MTP18 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MTP18 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse heart tissue. Positive IP detected in mouse heart tissue. Observed molecular weight by Western-blot: 16-18 kDa

Preparation

This antibody was obtained by immunization of MTP18 recombinant protein (Accession Number: NM_016498). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112793(MTP18 antibody) at dilution of 1:300



IP Result of anti-MTP18 (IP:Catalog No:112793, 4ug; Detection:Catalog No:112793 1:500) with mouse heart tissue lysate 4000ug.