

MYH10 antibody

Catalog Number: 112929

Product name

MYH10 antibody

IP: 1:500-1:5000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MYH10 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 200 kDa

Preparation

This antibody was obtained by immunization of MYH10 recombinant protein (Accession Number: NM_001256095). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112929(MYH10 antibody) at dilution of 1:1500



Immunofluorescent analysis of HepG2 cells, using MYH10 antibody Catalog No:112929 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



MYH10 antibody



Catalog Number: 112929

IP Result of anti-MYH10 (IP:Catalog No:112929,

4ug; Detection:Catalog No:112929 1:1500) with mouse brain tissue lysate 3440ug.