

NAT5 antibody

Catalog Number: 112963

Product name

NAT5 antibody

Immunogen

[Human NAT5 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NAT5 Rabbit Polyclonal antibody. Positive IP detected in HL-60 cells. Positive WB detected in HeLa cells, HL-60 cells, human brain tissue. Observed molecular weight by Western-blot: 18-20 kDa

Preparation

This antibody was obtained by immunization of NAT5 recombinant protein (Accession Number: NM_016100). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112963(NAT5 antibody) at dilution of 1:500



IP Result of anti-NAT5 (IP:Catalog No:112963, 4ug; Detection:Catalog No:112963 1:500) with HL-60 cells lysate 3040ug.