

# NBAS antibody

Catalog Number: 112965

## General Information

### Product name:

NBAS antibody

### Immunogen:

### Applications:

ELISA, WB, IF, IP

### Specificity:

Human, Mouse, Rat; other species not tested.

### Antibody description:

NBAS Rabbit Polyclonal antibody. Positive IP detected in SH-SY5Y cells. Positive WB detected in SH-SY5Y cells. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 220 kDa

## Properties

### Preparation:

This antibody was obtained by immunization of NBAS recombinant protein (Accession Number: BC051792). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal

**Ig Type:** Rabbit IgG

## Applications


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IF: 1:20-1:200

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:112965(NBAS antibody) at dilution of 1:500




IP Result of anti-NBAS (IP:Catalog No:112965, 4ug; Detection:Catalog No:112965 1:500) with SH-SY5Y cells lysate 1600ug.

# NBAS antibody



Catalog Number: 112965

---



Immunofluorescent analysis of (-20oc Ethanol) fixed SH-SY5Y cells using Catalog No:112965(NBAS Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)