

# Myocilin antibody



Catalog Number: 112989

## Product name

Myocilin antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Myocilin Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human stomach tissue, mouse skeletal muscle tissue. Positive IP detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 66 kDa

## Preparation

This antibody was obtained by immunization of My<sup>o</sup>Cilin recombinant protein (Accession Number: NM\_000261). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112989(MYoc antibody) at dilution of 1:300



IP Result of anti-MYoc (IP:Catalog No:112989, 4ug; Detection:Catalog No:112989 1:500) with mouse skeletal muscle tissue lysate 4000ug.