

NECAP2 antibody



Catalog Number: 113093

Product name

NECAP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NECAP2 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, human liver tissue, mouse cerebellum tissue. Observed molecular weight by Western-blot: 35-37 kDa

Preparation

This antibody was obtained by immunization of NECAP2 recombinant protein (Accession Number: NM_018090). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113093(NECAP2 antibody) at dilution of 1:600