

# NPY antibody

Catalog Number: 113226

---

## Product name

NPY antibody

## Immunogen

[Human NPY Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NPY Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer.

## Preparation

This antibody was obtained by immunization of NPY recombinant protein (Accession Number: NM\_000905). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113226(NPY antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113226(NPY antibody) at dilution of 1:50 (under 40x lens)