

NT5C3 antibody

Catalog Number: 113393

General Information

Product name:

NT5C3 antibody

Immunogen:

Applications:

ELISA, WB, IHC

Specificity:

human, mouse, rat; other species not tested.

Antibody description:

NT5C3 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer. Positive WB detected in human skeletal muscle tissue, mouse kidney tissue, mouse placenta tissue. Observed molecular weight by Western-blot: 33-38 kDa

Properties

Preparation:

This antibody was obtained by immunization of NT5C3 recombinant protein (Accession Number: BC015856). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal

Ig Type:Rabbit IgG


Applications

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113393(NT5C3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113393(NT5C3 antibody) at dilution of 1:50 (under 10x lens)