

Nucleostemin antibody



Catalog Number: 113410

Product name

Nucleostemin antibody

Specificity

Human; other species not tested.

Antibody description

Nucleostemin Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in K-562 cells, HEK-293 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 60-62 kDa

Preparation

This antibody was obtained by immunization of Nucleostemin recombinant protein (Accession Number: NM_014366). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113410(GNL3 antibody) at dilution of 1:1500



IP Result of anti-GNL3 (IP:Catalog No:113410, 4ug; Detection:Catalog No:113410 1:600) with HeLa cells lysate 4000ug.