

# OAS1 antibody



Catalog Number: 113458

## Product name

OAS1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

OAS1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of OAS1 recombinant protein (Accession Number: NM\_002534). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113458(OAS1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:113458(OAS1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:113458(OAS1

# OAS1 antibody



Catalog Number: 113458

---

antibody) at dilution of 1:100 (under 40x lens)