

OTC antibody



Catalog Number: 113510

Product name

OTC antibody

Immunogen

[Human OTC Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

OTC Rabbit Polyclonal antibody. Positive WB detected in rat liver tissue, mouse liver tissue. Observed molecular weight by Western-blot: 36-40 kDa

Preparation

This antibody was obtained by immunization of OTC recombinant protein (Accession Number: NM_000531). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



rat liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113510(OTC Antibody) at dilution of 1:1000