

# p115, USO1 antibody



Catalog Number: 113532

## Product name

p115, USO1 antibody

## Immunogen

[Human p115, USO1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

p115, USO1 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IHC detected in human gliomas tissue. Positive IP detected in mouse brain tissue. Positive WB detected in HepG2 cells, HEK-293 cells, human brain tissue, mouse testis tissue, mouse thymus tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 115kd

## Preparation

This antibody was obtained by immunization of p115, USO1 recombinant protein (Accession Number: NM\_003715). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IF, IP, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:20-1:200


## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113532(p115, USO1 Antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113532(p115, USO1 antibody) at dilution of 1:100 (under 10x




lens)

Immunofluorescent analysis of Hela cells, using USO1 antibody Catalog No:113532 at 1:50

# p115, USO1 antibody

Catalog Number: 113532


dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



DNA dye).

IP Result of anti-p115, USO1 (IP:Catalog No:113532, 3ug; Detection:Catalog No:113532

1:1500) with mouse brain tissue lysate 7000ug.



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:113532(p115, USO1 antibody) at dilution of 1:100 (under 40x lens)