

PAOX antibody



Catalog Number: 113574

Product name

PAOX antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PAOX Rabbit Polyclonal antibody. Positive WB detected in human stomach tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 56kd

Preparation

This antibody was obtained by immunization of PAOX recombinant protein (Accession Number: NM_152911). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:113574(PAOX antibody) at dilution of 1:500