

# PDRG1 antibody

Catalog Number: 113636

## Product name

PDRG1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PDRG1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IF detected in HepG2 cells. Positive IHC detected in human kidney tissue, human liver tissue. Observed molecular weight by Western-blot: 16 kDa

## Preparation

This antibody was obtained by immunization of PDRG1 recombinant protein (Accession Number: NM\_030815). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113636(PDRG1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113636(PDRG1 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human kidney using Catalog No:113636(PDRG1 antibody) at dilution of 1:100 (under 10x lens)

# PDRG1 antibody

Catalog Number: 113636

---

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using PDRG1 antibody Catalog No:113636 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).