

PCBD1 antibody

Catalog Number: 113683

Product name

PCBD1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

human, mouse, rat; other species not tested.

Antibody description

PCBD1 Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue. Positive IHC detected in human ovary tumor. Positive IF detected in Hela. Observed molecular weight by Western-blot: 12kd

Preparation

This antibody was obtained by immunization of PCBD1 recombinant protein (Accession Number: NM_000281). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:113683(PCBD1 antibody) at dilution of 1:100



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:113683(PCBD1 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded

PCBD1 antibody

Catalog Number: 113683

human ovary tumor using Catalog
No:113683(PCBD1 antibody) at dilution of 1:50
(under 40x lens)



Immunofluorescent analysis of HeLa cells, using
PCBD1 antibody Catalog No:113683 at 1:25
dilution and Rhodamine-labeled goat anti-rabbit
IgG (red).