

# PCDH11X antibody



Catalog Number: 113690

## Product name

PCDH11X antibody

## Specificity

Human; other species not tested.

## Antibody description

PCDH11X Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human brain tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 135 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_032968). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000


IHC: 1:20-1:200

IF: 1:10-1:100


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113690(PCDH11X antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human brain using Catalog No:113690(PCDH11X antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human brain using Catalog No:113690(PCDH11X antibody) at dilution of 1:50 (under 40x lens)

# PCDH11X antibody



Catalog Number: 113690

---



Immunofluorescent analysis of HepG2 cells, using PCDH11X antibody Catalog No:113690 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).