

PCYT2 antibody

Catalog Number: 113760

Product name

PCYT2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PCYT2 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, HepG2 cells, MCF7 cells, mouse testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of PCYT2 recombinant protein (Accession Number: NM_002861). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113760(PCYT2 antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:113760(PCYT2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)