

PDAP1 antibody



Catalog Number: 113762

Product name

PDAP1 antibody

IF: 1:10-1:100

Specificity

human, mouse, rat; other species not tested.

Antibody description

PDAP1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells. Positive IHC detected in human breast cancer. Positive IF detected in Hela. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of PDAP1 recombinant protein (Accession Number: NM_014891). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113762(PDAP1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113762(PDAP1 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human breast cancer using Catalog

PDAP1 antibody



Catalog Number: 113762

No:113762(PDAP1 antibody) at dilution of 1:50
(under 40x lens)



Immunofluorescent analysis of HeLa cells, using
PDAP1 antibody Catalog No:113762 at 1:25
dilution and Rhodamine-labeled goat anti-rabbit
IgG (red).