

# PHF11 antibody



Catalog Number: 113806

## Product name

PHF11 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PHF11 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse thymus tissue. Positive IHC detected in human lymphoma tissue. Observed molecular weight by Western-blot: 34kd

## Preparation

This antibody was obtained by immunization of PHF11 recombinant protein (Accession Number: NM\_001040444). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113806(PHF11 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:113806(PHF11 antibody) at dilution of 1:100 (under 10x lens)